

## Holding and Positioning Unit HPU

Holding and Positioning Unit for positioning FlexMEAs or NeuroNexus Probes connected to Miniature Preamplifiers  $\mu$ PA or MPA via Rod Mount



### Technical Specification

#### HPU

The holding and positioning unit HPU allows to fix the rod mount NN-RM in a stable position. A miniature preamplifier  $\mu$ PA or MPA is mounted to the rod mount, the preamplifier is connected to a FlexMEA or a NeuroNexus probe.

Note: The color of the HPU is variable!



**Important: The rod mount NN-RM and the miniature preamplifier need to be ordered separately.**

Type	Height with Base (mm)	Weight (g)	Diameter of the base (mm)
HPU160-S	234	1430	60
HPU160-L	198	1480	80
HPU290-S	364	1450	60
HPU290-L	328	1500	80